Version 01 – 01.26

EMPOWERING TOMORROW

en.siger.cn

Zhuhai Siger Medical Equipment Co., Ltd

Add : Building 2, No.1 Chuangxin Yi Road, Tangjiawan Town, Zhuhai City, Guangdong Province, P. R. China

Tel : +86 756 3881087 Email : contactus@siger.cn Web : en.siger.cn


: Siger medical : Siger_medical

SIGER

V1000

Leading Comfort Forward

V Series Dental Unit



Empowered by Expertise
Comfort Evolved

A Trusted Choice – Siger V1000 dental unit, brings superior performance for professional care, and human-centered design for enhanced comfort.

Boost Efficiency
Command with Ease





Whip Arm Dentist Element



Lightweight Pull Rod

V1000 whip arm version dentist element combined with lightweight pull rods, is easy to operate. Smooth in motion and stable in return, the ergonomic structure reduces hand and wrist fatigue for precise and efficient treatment.



Classic Control Panel

V1000 doctor control panel has clear button icons and excellent fingertip tactile feeling. The unit always responds promptly to your instructions.






Hanging-Type Dentist Element

Spacious instrument table with a newly upgraded control system and 5 instrument slots, allowing effortless operation and seamless transitions throughout the procedure.

Redefining Ergonomics & Comfort

Our dual design balances patient comfort with effortless operation for the dentist.




Comfortable Butterfly Backrest:

The wing-shaped arc naturally conforms to the body, providing comfortable and secure support for the back and arms.



Proprietary Foam Upholstery:

Cloud-like comfort is achieved through our precision-engineered foam density, thickness and shape.



Double-jointed Headrest:

Meet the needs of different treatment positions



SmoothControl™ Chair Comforting System:

Chair starts and stops softly, ensures smooth and gentle movement, providing next-level care to patients.




Low Minimum Chair Position:

Featuring a super low position, convenient to doctor and patient for safe and precise treatment.



Reinforced Chair Frame and Base:

Solid design eliminates chair sway during movements.



V1000 Upgraded with Extra-size Doctor Table

Large Tray, Six slots, Three Doctor Programs



User-Friendly Design Professional Assurance



Detachable Cuspidor & Water Outlet

Designed to eliminate hygiene blind spots.



Removable Suction Tubes & Filter

Enables efficient daily cleaning and maintenance



2L Large Bottle

Fewer refills, less interruption

Top-Filling Design

Refill in one step, simply and swiftly



V2 Pro Dental Light

Provides high-quality, virtually shadowless illumination for optimal visibility.




Light Source: LED
Color Temperature: 2700 - 5500 K
Illuminance: 8000 - 70000 LUX
Light Spot Size: 90*180mm (@700 MM)
Voltage: 12 - 24 Vac
Power: 18 VA

V1000 Assistant Element

Work in Sync: Enhancing Collaboration



Multifunction Foot Control



Let the Color Speak

*Colors shown are for reference only. For actual colors, please refer to the color swatches.

Standard Seamless Upholstery (PU Leather)

H1 Light Blue	H2 Jewelry Blue	H3 Grass Green	H4 Saffron Yellow	H5 Orange Red	H6 Light Grey	H7 Lake Green
H8 Pearl White	H9 Medium Blue	H10 Apple Green	H11 Bright Silver	H12 Purplish Red	H13 Chinese Red	H16 Black

Optional Sewn Upholstery (Silicone-coated Leather)

S1 Tundra	S2 Ciel	S3 Mocha	S4 Slate	S5 Truffle	S6 Terre	S7 Orange	S8 Bordeaux
S9 Ocean	S10 Raisin	S11 Black	S12 Fuscia	S13 Lilac	S14 Ivory	S15 Yellow	S16 Dark blue

Optional Italian Seamless Upholstery

T5003 Deep Blue	T1009D Tech Blue	T3491 Light Blue	T15464 Sky Blue	T3493 Lake Green	T15540 Emerald	T12143 Apple Green	T14766 Pearl Red	T7791 Orange
T6141 Bright Orange	T23337 Light Khaki	T3494 Light Gray	T23989 Silver	T3518 Silver Grey	T12504 Platinum	T7502 Dark Grey	T3913 Jet Black	