Version 01 – 01, 26

EMPOWERING TOMORROW

en.siger.cn

Zhuhai Siger Medical Equipment Co., Ltd

Add : Building 2, No.1 Chuangxin Yi Road, Tangjiawan Town, Zhuhai City, Guangdong Province, P. R. China


Tel : +86 756 3881087 Email : contactus@siger.cn Web : en.siger.cn

: Siger medical : Siger_medical

SIGER

U200
U100

U Series Dental Unit



A Trusted Partner
for Every Treatment



U200


A TRUSTED PARTNER FOR EVERY TREATMENT

Led by intelligent technology and built on exceptional quality, the U200 dental unit delivers a clinical experience that is comfortable, stable, and highly efficient.

U200


SMART CONTROL FOR PRECISION DENTISTRY

Powered by SigerSmart OS™, clinical workflows are reimagined as smart interactions—enabling dentists to deliver precise treatments with ease.



Real-Time Status Display

Delivers precise control over motor speed, forward/reverse rotation, and scaler power.




Personalized Configuration

Customizable parameters for intuitive and effortless operation.



Smart Synchronization

Customizable linking of cup rinse, cuspidor flush, dental light, and chair position for smooth, efficient operation.



Disinfection System

One-touch automatic pipeline disinfection. (Optional)

EVERY DETAIL THOUGHTFULLY COMFORTABLE

Where Patient Experiences a New Standard of Comfort.



Low Minimum Chair Position

Featuring a super low position, convenient to doctor and patient for safe and precise treatment.



Double-jointed Headrest

Meet the needs of different treatment positions

Comfortable Butterfly Backrest

The wing-shaped arc naturally conforms to the body, providing comfortable and secure support for the back and arms.





SmoothControl™ Chair Comforting System

Chair starts and stops softly, ensures smooth and gentle movement, providing next-level care to patients.

Reinforced Chair Frame and Base

Solid design eliminates chair sway during movements.

BIOTUBE™ AUTOMATIC PIPELINE DISINFECTION SYSTEM

Biotube™ automatic disinfection system delivers one-touch tube flushing and disinfection, effectively minimizing healthcare-associated infection risks.





- **Fully Visible Process**
  Disinfection progress is entirely transparent, giving operators clear control and confidence.

- **Automated Record Storage**
  Supports over 300 traceable records for rigorous and efficient management.



Fully automatic pipeline disinfection, professionally protecting every treatment.

EXQUISITE DESIGN

Take full control of every treatment.



Precise Temperature Control

Consistent temperatures for patient comfort.



Top-Filling Design

Refill in one step, simply and swiftly.

Detachable Cuspidor & Water Outlet

Designed to eliminate hygiene blind spots.



Removable Suction Tubes & Filter

Enables efficient daily cleaning and maintenance.



DESIGNED FOR EFFICIENT CLINICAL WORKFLOW

V3 Dental Light

Light Source: LED
Color Temperature: 2700 – 5500 K
Illuminance: 7500 – 90000 LUX
Light Spot Size: 90 × 180 mm (@700 mm)
Voltage: 12 – 24 V (AC/DC)
Power: 22 VA





V2 Pro Dental Light

Light Source: LED
Color Temperature: 3000 – 5500 K
Illuminance: 8000 – 70000 LUX
Light Spot Size: 90 × 180 mm (@700 mm)
Voltage: 12 – 24 V (AC/DC)
Power: 18 VA

Optional Backrest Designs

designed to suit different clinical preferences



Butterfly-shaped Backrest ▶
Wider design for enhanced support and comfort.

◀ Fan-shaped Backrest
Streamlined design for greater simplicity and freedom of movement.



LET THE COLOR SPEAK

*Colors shown are for reference only. For actual colors, please refer to the color swatches.

Standard Seamless Upholstery (PU Leather)

H1 Light Blue	H2 Jewelry Blue	H3 Grass Green	H4 Saffron Yellow	H5 Orange Red	H6 Light Grey	H7 Lake Green
H8 Pearl White	H9 Medium Blue	H10 Apple Green	H11 Bright Silver	H12 Purplish Red	H13 Chinese Red	H14 Black

Optional Sewn Upholsterys (silicone-coated leather)

S1 Tundra	S2 Ciel	S3 Mocha	S4 Slate	S5 Truffle	S6 Terre	S7 Orange	S8 Bordeaux
S9 Ocean	S10 Raisin	S11 Black	S12 Fuscia	S13 Lilac	S14 Ivory	S15 Yellow	S16 Dark blue

Optional Italian Seamless Upholstery

T5003 Deep Blue	T1009D Tech Blue	T3491 Light Blue	T15464 Sky Blue	T3493 Lake Green	T15540 Emerald	T12143 Apple Green	T14766 Pearl Red	T7791 Orange
T6141 Bright Orange	T23337 Light Khaki	T3494 Light Gray	T23989 Silver	T3518 Silver Grey	T12504 Platinum	T7502 Dark Grey	T3913 Jet Black	